


Quick Start Guide: Lymph2Cx Assay DLBCL cell-of-origin Classifier Web Site

Home page: <https://llmpp.nih.gov/LYMPHCX/index.shtml>



LYMPHCX Home

<https://llmpp.nih.gov/LYMPHCX/>

Lymphoma/Leukemia
Molecular Profiling Project
(LLMPP)

**Lymph2Cx Assay
DLBCL cell-of-origin Classifier**

Page Updated: Wednesday, 20-Aug-2014 08:22:55 EDT

Access to the DLBCL subtype Classification tool is restricted to a limited set of users.

- [Access the Classification Tool using your NIH/Federated login](#)

or

- [Access the Classification Tool using a local login/password](#)

Disclaimer: This test is designed to be used on FFPE samples that have been confirmed by Pathologists to be of the DLBCL subtype with at least 60% tumor content, and have been processed according to the Nanostring SOP. Other samples may produce unexpected or inaccurate results. Although some quality control is incorporated into the algorithm it is recommended that users additionally exclude samples that have been flagged by the NanoString nSolver Analysis software as being of poor quality, or otherwise appear as outliers. This tool has not been evaluated by the FDA or any other regulatory agency, and is not intended for clinical use.

Direct questions/support issues to llmpp-support@bimas.cit.nih.gov

 NIH Bioinformatics support provided by [BIMAS/CBEL/CIT](#).
We can be contacted by [email](mailto:llmpp-support@bimas.cit.nih.gov).




At the bottom of the every page you will find an email link to llmpp-support@bimas.cit.nih.gov . Use this email for any questions, problems or requests related to this Web Site.

- Accessing the Classification Tool with an NIH Login and Password


Users who have an NIH Login and Password should click the “Access the Classification Tool using your NIH/Federated login” link.

- [Access the Classification Tool using your NIH/Federated login](#)
- or
- [Access the Classification Tool using a local login/password](#)

You will then be directed to the NIH Single Sign On site to authenticate. For the Account Type value, select “NIH Staff”.



You may authenticate using either your NIH User Name/Password or with your PIV card/PIN. To authenticate with you NIH User Name and Password, enter your NIH User Name and Password and click the “Log in” button.



- Accessing the Classification Tool with a provided login and password

Users without an NIH Login and Password have been provided with a user name and password and should click on the “Access the Classification Tool using a local login/password” link.

- [Access the Classification Tool using your NIH/Federated login](#)
- or
- [Access the Classification Tool using a local login/password](#)

Depending on your browser, you will be prompted to “Log in” or “Authenticate”. The example below is using the Chrome Browser. The first time you access the site, enter the user name and password information that was provided to you and click the “Log In” button.

Authentication Required

The server https://lmpp.nih.gov:443 requires a username and password. The server says: Enter your LYMPHCX User Name and Password.

User Name:

Password:

You will then be required to change the password

You must reset your Lymph2Cx Web Site password.

A good password should have a mix of upper/lower case letters, numbers and special characters such as @ # 0 * ? . The Lymph2Cx password must be between 5 and 8 characters long.

Reset Local Password for Iso_test

New Password: (5 to 8 characters long)

Confirm New Password:

NOTE: On successful resetting of your local password, you will be immediately prompted to "log in" or "Authenticate" with your user name and new password.

After submitting your new password, you will be required to “Log in” or “Authenticate” again with your user name and new password.

- Uploading Lymph2Cx Assay data

Lymph2Cx Assay DLBCL cell-of-origin Classifier

User John Powell as jpowell


Input data should be extracted from nSolver and saved as a tab delimited text file. It should include columns for Gene Name, Accession #, and Class name, remaining columns should be unnormalized, untransformed data headed by the sample identifier.

An example input file in the correct format can be found [here](#).

Select your File and Upload for Classification

Data file: No file chosen

Click on “Choose File” and navigate to and select the TAB delimited file extracted from NanoString’s nSolver. In this example, “input_example.tab” has been selected. Click the “Open” button to complete the selection.



Click the “Upload” button to upload the selected file to the Lymph2Cx assay cell-of-origin classifier

Select your File and Upload for Classification

Data file: input_example.tab

- Lymph2cx Assay cell-of-origin Classifier Results Display

Formatted for

← Form for Retrieving Results

Order	Samples	Normalization	ABC.likelihood	Prediction
1	Sample.1	190	0.01	GCB
2	Sample.2	1831	0.03	GCB
3	Sample.3	438	1	ABC
4	Sample.4	79	0	GCB
5	Sample.5	307	0	GCB
6	Sample.6	45	0	GCB
7	Sample.7	247	1	ABC
8	Sample.8	144	0	GCB
9	Sample.9	210	0.19	intermediate/Unclassified
10	Sample.10	276	0.32	intermediate/Unclassified
11	Sample.11	129	0.93	ABC
12	Sample.12	64	0.84	intermediate/Unclassified

← Display of Results

Column definition:

← Column Definitions

Samples

Sample identifier extracted from the header of the input file

Use the box at the top of the display to retrieve the prediction results formatted as TAB delimited text, for PDF or for Excel.

Formatted for
Excel
PDF
Tab delimited text

The lines in the TAB delimited retrieved files are terminated with a "linefeed" character. If you use Notepad to view the TAB delimited text, the lines will display "run together". Use WordPad or Word to view the files.

- Lymph2cx Assay cell-of-origin Classifier exported PDF example

Classification results for: input_example.tab

**Lymph2Cx Assay
DLBCL cell-of-origin Classifier**

John Powell as nih.gov:jip

Portal Code Rev 258 2014-09-03 10:09:26 -0400 (Wed, 03 Sep 2014)
lymphcx.R Rev 200 2014-05-14 12:48:57 -0400 (Wed, 14 May 2014)
param.txt Rev 189 2014-05-09 10:10:08 -0400 (Fri, 09 May 2014)
gene.id Rev 189 2014-05-09 10:10:08 -0400 (Fri, 09 May 2014)

Order	Samples	Normalization	ABC Likelihood	Prediction
1	Sample.1	190	0.01	GCB
2	Sample.2	1831	0.03	GCB
3	Sample.3	438	1	ABC
4	Sample.4	79	0	GCB
5	Sample.5	307	0	GCB
6	Sample.6	45	0	GCB
7	Sample.7	247	1	ABC
8	Sample.8	144	0	GCB
9	Sample.9	210	0.19	intermediate/Unclassified
10	Sample.10	276	0.32	intermediate/Unclassified
11	Sample.11	129	0.93	ABC
12	Sample.12	64	0.84	intermediate/Unclassified

Column definition:

Samples

Sample identifier extracted from the header of the input file

Normalization

Summary of housekeeping genes. Higher values generally indicate improved quality. Samples with a normalization of less than 20 are deemed to be poor quality

ABC Likelihood

Estimated probability that the sample is ABC

Prediction

Final call for sample

Disclaimer: This test is designed to be used on FFPE samples that have been confirmed by Pathologists to be of the DLBCL subtype with at least 60% tumor content, and have been processed according to the Nanostring SOP. Other samples may produce unexpected or inaccurate results. Although some quality control is incorporated into the algorithm it is recommended that users additionally exclude samples that have been flagged by the NanoString nSolver Analysis software as being of poor quality, or otherwise appear as outliers. This tool has not been evaluated by the FDA or any other regulatory agency, and is not intended for clinical use.

Direct questions/support issues to lmpp-support@bimas.cit.nih.gov

Appendix A

- What ifs

What if I forget my local user name or password?

You will need to send an email to lmpp-support@bimas.cit.nih.gov . Reset password information will be send to the email of record when the account was created.

What if I want to change my password?

You will need to send an email to lmpp-support@bimas.cit.nih.gov . Reset password information will be send to the email of record when the account was created.

What if I want to log out of the Lymph2Cx Portal?

If you logged in with an NIH/Federated login, you can click the “Log out” link near the top of the page. If you logged in with a local login, the only way to log out is to completely close your browser.

What if my uploaded Lymph2Cx Assay data file does not meet the requirements?

The Lymph2Cx classified system does a thorough job of checking the input file for any issues before attempting to run the classification. Any issues detected with the file are reported back to you.

What if the classifier cannot determine the cell-of-origin (COO)?

If the input was deemed to be to poor quality to attempt classification, the classifier prediction would be “Poor Quality”. If quality is Okay, the classifier returns “Intermediate/Unclassified” for those samples cannot be classified ABC or GCB.

Appendix B

- Required format of the Lymph2Cx Assay file extracted from nSolver Analysis Software

- The file must be a TAB-separated text file
- The first line of the file must be a header line. The first 3 header columns must be 'Gene Name', 'Accession #', and 'Class Name'. The remaining header columns are sample identifiers.
- The nSolver measurement values are expected to be zero or positive numbers. Negative values and missing measurements are not allowed.
- The following list of twenty Gene Name & Accession values must be present in the file:

Gene Name	Accession #
ASB13	NM_024701.3
CCDC50	NM_174908.3
CREB3L2	NM_194071.2
CYB5R2	NM_016229.3
IRF4	NM_002460.1
ISY1	NM_020701.2
ITPKB	NM_002221.3
LIMD1	NM_014240.2
MAML3	NM_018717.4
MME	NM_000902.2
MYBL1	XM_034274.14
PIM2	NM_006875.2
R3HDM1	NM_015361.2
RAB7L1	NM_001135664.1
S1PR2	NM_004230.2
SERPINA9	NM_001042518.1
TNFRSF13B	NM_012452.2
TRIM56	NM_030961.1
UBXN4	NM_014607.3
WDR55	NM_017706.4

- Example Lymph2Cx Assay File containing 6 samples

Gene Name	Accession #	Class Name	Sample.1	Sample.2	Sample.3	Sample.4	Sample.5	Sample.6
ASB13	NM_024701.3	Endogenous	212	3639	111	876	617	42
CCDC50	NM_174908.3	Endogenous	634	4898	1360	242	345	205
CREB3L2	NM_194071.2	Endogenous	29	1401	805	19	91	180
CYB5R2	NM_016229.3	Endogenous	150	2236	34	189	37	8
IRF4	NM_002460.1	Endogenous	1349	11641	1934	291	794	831
ISY1	NM_020701.2	Endogenous	121	1280	208	71	230	22
ITPKB	NM_002221.3	Endogenous	823	6800	690	227	1623	37
LIMD1	NM_014240.2	Endogenous	221	1275	403	22	105	112
MAML3	NM_018717.4	Endogenous	65	194	21	151	75	45
MME	NM_000902.2	Endogenous	390	2699	5	99	36	536
MYBL1	XM_034274.14	Endogenous	190	3547	19	101	54	15
PIM2	NM_006875.2	Endogenous	7140	23675	2191	483	1228	3721
R3HDM1	NM_015361.2	Endogenous	310	2558	578	99	449	63
RAB7L1	NM_001135664.1	Endogenous	493	5748	4646	319	385	163
S1PR2	NM_004230.2	Endogenous	168	2059	144	255	540	341
SERPINA9	NM_001042518.1	Endogenous	326	998	1	330	76	1528
TNFRSF13B	NM_012452.2	Endogenous	150	2653	1490	21	192	1
TRIM56	NM_030961.1	Endogenous	49	767	280	16	110	9
UBXN4	NM_014607.3	Endogenous	898	4412	986	208	823	294
WDR55	NM_017706.4	Endogenous	150	1857	484	130	290	52
NEG_A	ERCC_00096.1	Negative	0	0	0	0	1	1
NEG_B	ERCC_00041.1	Negative	0	1	1	1	2	0
NEG_C	ERCC_00019.1	Negative	1	4	1	1	0	0
NEG_D	ERCC_00076.1	Negative	0	7	1	1	0	1
NEG_E	ERCC_00098.1	Negative	1	0	2	2	1	2
NEG_F	ERCC_00126.1	Negative	0	2	0	0	0	0
POS_A	ERCC_00117.1	Positive	26180	30754	31312	34199	30995	20442
POS_B	ERCC_00112.1	Positive	6221	7901	7041	7836	7998	5102
POS_C	ERCC_00002.1	Positive	1538	1708	2155	2320	2068	1257
POS_D	ERCC_00092.1	Positive	535	536	558	498	572	325
POS_E	ERCC_00035.1	Positive	99	192	144	165	138	81
POS_F	ERCC_00034.1	Positive	22	34	35	31	28	18